



# Fit With Clean Shield ZERO™ Less Touch, More Protection

Reducing Cross-Contamination While Improving  
Glove Dispensing Efficiency



Fit  
With Clean Shield ZERO™

**PRIMED'S CLEAN SHIELD ZERO™** is designed to help clinicians dispense gloves quickly and hygienically, eliminating contact with the film, reducing cross-contamination risks. Traditional glove boxes often require touching the film, increasing contamination risk, while small openings make access difficult, especially as the box empties.

With a wider, rotated opening, Clean Shield ZERO™ ensures efficient, one-at-a-time dispensing, improving workflow and accessibility until the last glove. Its design eliminates unnecessary contact, wasted gloves, and hassle, enhancing safety and ease of use.



## KEY BENEFITS FOR CLINICIANS

- **Enhanced Infection Control**  
Reduced risk of cross-contamination.
- **Increased Efficiency**  
Quick, easy dispensing with minimal disruption.
- **Less Waste**  
Gloves dispense one at a time, preventing unnecessary usage.
- **Improved Accessibility**  
Easy access to gloves, even when the box is near empty.

## Experience the Difference

With Clean Shield ZERO™, PRIMED continues to innovate for a safer and more effective healthcare environment.